



Professional Orthopedic Products

Presents the VISSION Stretch-UP

The Stretch-UP is a sock intended to be used when resting or at night time. The sock is designed to keep tension on the tissue (plantar fascia) so it heals in a stretched position at night.

When properly applied, the sock will hold the forefoot and ankle joints in a slight upward position, thereby preventing the plantar fascia from retracting and healing in a contracted position. In the morning the plantar fascia is at its maximum stretch, which reduces or eliminates the sharp pain that often occurs during the first steps.


Bilateral design fits both left and right foot.

Indications:

- Plantar Fasciitis
- Heel Spur

<u>Item #</u>	<u>Size</u>	<u>Circumference*</u>
341080011	Small	< 20 cm
341080012	Medium	20 - 40 cm
341080013	Large	> 40 cm

* Calf circumference



Professional Orthopedic Products

Toll Free Phone: 1-866-334-4441

Toll Free Fax: 1-888-484-4431

Email: csr@poortho.com

www.poortho.com