

V-Strap™ Wrist Support



When you want to provide your patients with the immobilization they need and the easy application and comfort they demand in order to be compliant, prescribe the **V-Strap** Wrist Support.

Optimum Comfort

Made with perforated suede and a polypropylene felt liner, the V-Strap allows moisture to evaporate away from the skin. The anatomically contoured design and the adjustable thumb strap provide a comfortable fit and allow unrestricted motion of the fingers.

Superior Immobilization

The proximal and distal closure straps work in concert with the palmar and the dorsal aluminum stays to resist flexion and extension of the wrist. The blue palmar pad enhances this support by matching the contour of the depression in the palm.

V-Strap Retention System

Provides extraordinary ease of application since the "V" formed by the closure strap and the hook fastener prevents the strap from disengaging with the plastic buckle when applying or removing the wrist support.

Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Post operative care



8" wrist support



10 1/2" wrist & forearm support

Suggested HCPCS: L3908

Size	8" black left	8" black right	8" beige left	8" beige right
pediatric	223249	223259		
x-small	223241	223251	223141	223151
small	223242	223252	223142	223152
medium	223244	223254	223144	223154
large	223245	223255	223145	223155
x-large	223246	223256	223146	223156
xx-large	223247	223257		

Size	10 1/2" black left	10 1/2" black right	10 1/2" beige left	10 1/2" beige right
small	223262	223272	223162	223172
medium	223264	223274	223164	223174
large	223265	223275	223165	223175
x-large	223266	223276	223166	223176